


Jitron®

NON CONTACT FOREHEAD THERMOMETER

JTMI-603M


QUICK GUIDE



Thank you for choosing Jitron Non Contact Forehead Thermometer JTMI-603M. The thermometer does not radiate any infrared rays and is designed for non contact forehead temperature measurement. This Quick Guide provides important steps that help you to test your temperature. Please read the detailed User's Guide in the box carefully before using the meter.

Before you begin, make sure that the sensor is clean and without damage. Patients and thermometer should stay in steady state room condition for at least 30 minutes before the measurement. We recommend you measure 3 times and take the average reading as your temperature.

FOREHEAD MEASUREMENT



1. Forehead measurement is convenient. It requires a clean forehead to scan. Clean the forehead before the measurement.
2. Press the "SCAN" button to switch ON the meter. Aim the sensor at the forehead 0.5 to 3 cm away. Then press and hold the "SCAN" button until you hear a beep sound.
3. The measurement will be displayed in one second.
4. Wait 20 seconds for the next measurement.
5. As the forehead temperature is likely to be affected by many factors including sweat, oil, cosmetics and surroundings the forehead reading can be taken as a reference only.

MODE SELECTION

In power on mode, press the "MODE" button to switch between the modes HUMAN, OBJECT and both in silent modes. Set at Human mode for forehead temperature measurement.


°C / °F SWITCH

Degree Celsius (°C) or Fahrenheit (°F) unit display is available in this device. While the unit is ON, press and hold both the "MODE" button and the "MEMORY" button for approximately 3 seconds. This will change the mode to either °C or °F.

MEMORY


While in power ON mode, press & release "MEMORY" button to see the temperature stored in memory.


CLEANING THE LENS



Clean the sensor periodically to ensure an accurate reading. Please use a cotton swab with 70% concentrated alcohol to clean the lens. Allow the lens to dry for at least 1 minute before the next measurement or use a tissue paper to dry it.

BATTERY REPLACEMENT

The low battery icon  on the screen indicates low battery. Please open the battery cover by sliding the battery cover out. Replace with new AAA 1.5V alkaline batteries.



BATTERY:
AAA 1.5V ALKALINE

POWER

The Power will be OFF automatically if not in use for approximately 30 seconds.

ADDITIONAL FEATURES

Please refer to the User's Guide to understand object temperature measurement and silent mode.

Jitron®
S I N G A P O R E



Manufacturer:

Jitron Pte Ltd

8 Ubi Road 2, #05-01

Zervex, Singapore 408538

Tel: (65) 6748 6815

E-mail: enquiry@jitron.com

Web: www.jitron.com

Rev. 02/20