


Jitron®

NON CONTACT FOREHEAD THERMOMETER

JTMI-603M


QUICK GUIDE



Thank you for choosing Jitron Non Contact Forehead Thermometer JTMI-603M. The thermometer does not radiate any infrared rays and is designed for non contact forehead temperature measurement. This Quick Guide provides important steps that help you to test your temperature. Please read the detailed User's Guide in the box carefully before using the meter.

Before you begin, make sure that the sensor is clean and without damage. Patients and thermometer should stay in steady state room condition for at least 30 minutes before the measurement. We recommend you measure 3 times and take the average reading as your temperature.

FOREHEAD MEASUREMENT



1. Forehead measurement is convenient. It requires a clean forehead to scan. Clean the forehead before the measurement.
2. Press the "SCAN" button to switch ON the meter. Aim the sensor at the forehead 0.5 to 3 cm away. Then press and hold the "SCAN" button until you hear a beep sound.
3. The measurement will be displayed in one second.
4. Wait 20 seconds for the next measurement.
5. As the forehead temperature is likely to be affected by many factors including sweat, oil, cosmetics and surroundings the forehead reading can be taken as a reference only.

MODE SELECTION

In power on mode, press the "MODE" button to switch between the modes HUMAN, OBJECT and both in silent modes. Set at Human mode for forehead temperature measurement.

