



JBPM-903A


Digital Arm Blood Pressure Monitor




Regular to Large




2 Users
Start Key Design



Hypertension
Range Indication



90 Memory
Sets Each



Average
Memory



Irregular Heartbeat
Detection



Low Battery
Indicator



Alarm Clock
Setting



Date and Time
Settings





Arm Cuff
22~42cm



Large LCD Display
and Backlight

Irregular
Heartbeat Detection

Hypertension
Indication



Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that **the cuff is at the same level as your heart.** Relax your left hand.

Authorised Distributor: