



## Deluxe ARM TYPE

JBPM-903A



Regular to Large

- 2 Users x 90 Sets Memory
- Blood Pressure Range Indicator
- Irregular Pulse Detection

## Popular ARM TYPE

JBPM-902A



- 60 Sets Memory

Irregular Pulse Detection



Blood Pressure Range Indicator

Irregular Pulse Detection



## BLOOD PRESSURE MONITOR

### Premium ARM TYPE

BPI-901A



- 2 Users x 50 Sets Memory

### Traveler WRIST TYPE

BPI-801W  
Fuzzy Logic Technology

- 85 Sets Memory



## THERMOMETER

### NON CONTACT FOREHEAD THERMOMETER

JTMD-201M



- 1 Second measurement
- °C / °F selection
- Contactless

JTMI-603M

## DIGITAL THERMOMETER



- 10 Seconds measurement
- °C / °F selection
- Flexible tip

### EAR/FOREHEAD THERMOMETER



- 1 Second ear measurement
- No probe cover requirement
- °C / °F selection

JTMI-602M

## COMPRESSOR NEBULIZER

- Nebulizer Capacity: 8 ml
- Nebulisation Rates: 0.4 ml / minutes
- Sound Level: 52dBA
- Particle Size: 3.16 µm (MMAD)



JN-702C

## DiabeCHECK® BLOOD GLUCOSE MONITORING SYSTEM

DC-302MS

Test Strips:  
25's, 50's and 100's

- Tiny 1µl Blood Sample
- Result in < 6 Seconds
- Strip/Lancet Ejection Feature
- 300 Memory with Date & Time
- mg/dL or mmol/L (Region dependent)



## PROFESSIONAL STETHOSCOPE


- Adult Stainless Steel
- Adult Chrome Plated
- Pediatric Stainless Steel
- Cardiology Stainless Steel



- Tunable diaphragm
- Low and high frequency acoustic sensitivity


## LED PENLIGHT

JT-907



Penlight

CF-102RB



Cobalt Filter (Red & Blue)  
Accessory for Eye Testing

## DIGITAL (MERCURY-FREE) SPHYGMOMANOMETER

JT-E301A



- Mercury-Free Operation
- Manual Inflation and Air Release
- To be used with Stethoscope
- Portable

## OxyCHECK® FINGERTIP PULSE OXIMETER



JPO-707A

- SpO<sub>2</sub>
- Pulse Rate
- Perfusion Index

Manufacturer:

Jitron Pte Ltd

8 Ubi Road 2, #05-01  
Zervex, Singapore 408538  
Tel: (65) 6748 6815  
E-mail: enquiry@jitron.com  
Web: www.jitron.com

