



WS-CIC
(Complete In The Canal)



WS-ITC
(In The Canal)



WS-CIC & WS-ITC are ultra compact size with medium power output. These are discreet, yet still visible within the outer ear.

WS-20E is a compact size hearing aid, tiny shell design, medium power output & volume control.

WS-20E
(In The Ear)



WS-95H
(Behind The Ear)



WS-95E
(Behind The Ear)



WS-92EH
(Behind The Ear)



WS-95H, WS-95E, WS-92EH are Behind The Ear (BTE) type hearing aids with plastic housing design of different sizes, volume control & power on-off switch on selected models.

WS-3PE
(Pocket)



WS-3PE is a Pocket Type or Body-Aid hearing aid that can be carried in the pocket or put on the belt or waistband. It has high power output & volume control suitable for people with severe hearing loss and uses normal AA battery for power.



JITRON HEARING AID - SCALE



JITRON HEARING GUIDE

MODEL	SIZE	MILD	MODERATE	SEVERE	PROFOUND
1 WS-CIC - Complete In The Canal	Miniature	←→	←→		
2 WS-ITC - In The Canal	Tiny	←→	←→		
3 WS-20E - In The Ear	Mini	←→			
4 WS-95H - Behind The Ear	Small		←→	←→	
5 WS-95E - Behind The Ear	Small	←→	←→	←→	
6 WS-92EH - Behind The Ear	Medium	←→	←→	←→	
7 WS-3PE - Pocket Type	Pocket Size		←→	←→	

Authorised Distributor :