

Supreme



Available in seat widths :
40/46/51/56 cm
and choice of 3

different frame colours :

Blue/Maroon/Silver

Folded Width : 42 cm

Maximum Load : 136kg

Frame : Aluminium

Upholstery : Black Nylon

Armrest : Adjustable & Detachable

Front Castor : 200*50PU

Tyre N.W 63.kg

**ALL POWER WHEELCHAIRS
ARE FITTED WITH THE WHEELCHAIR
CONTROL SYSTEM FULLY DESIGNED AND
DEVELOPED IN SINGAPORE**

Extraordinary
Power Wheelchairs



Savers



Available in seat widths :
40 or 46 cm and
with standard frame colour : Silver/Blue

Folded Width : 42cm

Maximum Load : 113kg

Frame : Steel, Upholstery : Black Nylon

Armrest : Flip-Back & Non-Detachable

Front Castor : 200*50PU TYRE N.W 55Kg

Seat Width : 30/35 cm

Folded Width : 42 cm,

Frame : Aluminium, Color : Purple

Maximum Load : 50kg

Upholstery : Black Nylon

Armrest Height : Adjustable &

Detachable

Footrest : Swing-Away/Aluminium

Footplate : Plastic

Front Castor : 200*50PU Tyre N.W 36 kg

Junior



ACCESSORIES & REPLACEMENT PARTS



36AH BATTERY BOX 12793
17AH BATTERY BOX 12792



17AH 12790
36AH 12791



BATTERY CHARGER 12800



UNIVERSAL BATTERY TRAY 12801



PU TYRE ONLY 317 x 50mm



COMPLETE PU CASTOR 200 x 50mm



COMPLETE HUB WHEEL 12803



EXTERNAL MOTOR ASSEMBLY 12805



V1 CONTROL SYSTEM 12806

Authorised Distributor :