



STETHOSCOPE



SPECIFICATIONS

FEATURES

- Durable Tubing**
- Solid Chestpiece**
- Tunable Diaphragm**
- Nonchill Rim & Diaphragm**
- Double-leaf Binaural Spring inside the Tubing**
- Low and High Frequency Acoustic Sensitivity**

- Chestpiece: Dual Head Nonchill Rim & Diaphragm**
- Diaphragm: Tunable Diaphragm**
- Binaural: Dual-leaf Spring encased at 15° angle**
- Operating Temperature Range: 0°C ~ 40°C (32°F ~ 104°F)**
- Operating Humidity: Less than 95% RH**
- Storage Temperature Range: -10°C ~ 50°C (14°F ~ 122°F)**
- Storage Humidity: Less than 95% RH**

Authorised Distributor :

MODELS	TUBING	LENGTH	WEIGHT
Adult Stainless Steel JT-S601PF	Single-lumen	76cm	150g
Adult Chrome Plated JT-601PF	Single-lumen	76cm	130g
Cardiology Stainless Steel JT-S747PF	Double-lumen	69cm	215g
Pediatric Stainless Steel JT-S606PF	Single-lumen	76cm	125g