

# Jitron®

Biomedical & HealthCare  
SINGAPORE

www.jitron.com



# STETHOSCOPE



## FEATURES

- ★ Durable tubing
- ★ Solid chestpiece
- ★ Tunable diaphragm
- ★ Nonchill rim and diaphragm
- ★ Double-leaf binaural spring inside the tubing
- ★ Low and high frequency acoustic sensitivity

## SPECIFICATIONS

- Chestpiece : Dual head nonchill rim & diaphragm  
Diaphragm : Tunable diaphragm  
Binaural : Dual-leaf spring encased at 15° angle  
Operating temperature range : 0° C ~ 40° C  
(32°F ~ 104° F)  
Operating humidity : Less than 95% RH  
Storage temperature range : -10° C ~ 50° C  
(14°F ~ 122° F)  
Storage humidity : Less than 95% RH

Authorised Distributor :

MODELS	TUBING	LENGTH	WEIGHT
Adult Stainless Steel JT-S601PF	Single-lumen	76 cm	150 g
Adult Chrome Plated JT-601PF	Single-lumen	76 cm	130 g
Cardiology Stainless Steel JT-S747PF	Double-lumen	69 cm	215 g
Pediatric Stainless Steel JT-S606PF	Single-lumen	76 cm	125 g

Professional Stethoscopes for the healthcare industry from Jitron Singapore