

# Jitron®



[www.jitron.com](http://www.jitron.com)

*Biomedical & HealthCare, Singapore*  
**A Trusted Partner for Reliable Biomedical Products & Services**



## Jitron® Bio-Support Equipment Solutions

### Wheel Chairs

#### Specialised



Powerd & Joystick controlled hub motors with built-in electro-mech brakes.  
 Travel range of 28 km with fully charged 36AH Batteries

Seat width : 16", 18", 20"  
 Max. Load : 136 kg  
 Weight : Approx 63 kg

➤ **MOTORISED ALUMINUM - SUPREME**



Powerd & Joystick controlled external motors as standard with option to upgrade to hub motors

Seat width : 16" or 18"  
 Max. Load : 113 kg  
 Weight : Approx 55 kg

➤ **MOTORISED STEEL - SAVERS**



➤ **STAIR CLIMBER**

#### Aluminum



Standard Alum wheelchair with flip away armrests & detachable footrests. 24" PU mag wheels

Seat width : 16" or 18"  
 Weight : Approx 30lbs (14kg)

➤ **FLIP AWAY ARMRESTS & DETACHABLE FOOTREST - 4004/5**



Standard Alum wheelchair. Fixed armrest & footrests. Folding back with brake handle for carer. 24" PU mag wheels

Seat width : 16" or 18"  
 Seat to floor height : 20" (51 cm)  
 Weight : Approx 30lbs (14kg)

➤ **FIXED ARMRESTS & FOOTREST - 0708Mag**



Aluminum Push Chair with folding back. Fixed Armrest with removable footrests.

Seat width : 17" (43cm)  
 Seat to floor height : 18" (45cm)  
 Weight : Approx 22lbs (10kg)

➤ **ALUMINUM PUSH CHAIR-0707LT**

#### Steel, Chrome Plated



Removable armrests & swingaway removable footrests with rear 24" Solid Mag wheels.

Seat Width : 16" or 18"  
 Seat to Floor Height : 20" (51 cm)  
 Weight : Approx 44 lbs (20 kg)

➤ **DETACHABLE ARMRESTS/FOOTREST WHEELCHAIR - 0702/MAG**



Fixed armrests & Fixed footrest with 24" rear mag wheels. Economical and low in maintenance.

Seat Width : 18" (46 cm)  
 Seat to Floor Height : 20" (51 cm)  
 Weight : Approx 44 lbs (20 kg)

➤ **STANDARD WHEELCHAIR - 0700/MAG**



Standard Push Chair. Durable and long lasting in steel chrome. With fold down back for easy storage.

Seat width : 17" (43cm)  
 Seat to floor height : 20" (51 cm)  
 Weight : Approx 28lbs (13kg)

➤ **STANDARD PUSH CHAIR-0707FB**