





JBPM-902A DIGITAL BLOOD PRESSURE MONITOR

Arm Type



Arm
Measurement

Memory
60



Average
Memory Avg

Last **3** memory



Regular to
Large

Arm Cuff
22~42cm



Measuring Range :

Blood Pressure : 30~280 mmHg

Pulse Rate : 40~199 beats/min

Calibration Accuracy:

Blood Pressure : ± 3 mmHg

Pulse rate : $\pm 4\%$ of reading

Operating Environment :

10°C~40°C (50°F~104°F)

with relative humidity up to 85% (non condensing)

Storage/Transportation Environment :

-20°C~+50°C (-4°F~+122°F)

with relative humidity up to 85% (non condensing)

Battery : 4 x 1.5V AA alkaline batteries

Battery Life: Approx. 300 operations

Dimensions : 130mmX150mmX65mm (WXHXD)

Weight: 474g (including batteries and cuff)

Cuff Circumference : Approx. 22~42cm (Regular to Large)

LARGE DISPLAY

EASY TO USE

60 SETS MEMORY



Authorised Distributor :