



Portable Ultrasonic

JN-701U

NEBULIZER

Quiet

- Extremely quiet operation (noise level below 40dB)

Fast

- Fast flow-rate of up to 0.7ml/min

Safe

- Overheat protection
- Low battery indication

Effective

- Adjustable flow volume
- Particulate size of 0.5 to 5 microns increases medication absorption into the lungs

Lightweight

- Weight only 198g



Quiet

Fast

Safe

Effective

Lightweight



Portable Design

Accessories



Carrying Bag



Adult / Pediatric Masks



Disposable Ampoules



Spare Air Filters



AC/DC Adapter



Mouth Piece



Car Connection Cable



Rechargeable Battery

Specifications

Power	AC / DC Adapter (230VAC, 50Hz; 15VDC, 1200mA)
Nebulizer Power	12VDC, 800mA
Power Consumption	Approx. 12W
Sound Level	Below 40dB
Particle Range(MMAD)	0.5~5.0µm
Ampoule Capacity	Max. 5ml
Flow Rate	Max. 0.7ml / minute
Intermittent Use	20 min ON / 45 min OFF
Frequency	2.8MHz
Operation Conditions	10°C ~40°C (50°F~104°F); RH 10~95%
Storage Conditions	-25°C ~ 60°C (-13°F ~ 140°F); RH 10~95%
Dimensions	Approx. 85 mm x 55 mm x 130 mm (w/o batt.) Approx. 85 mm x 55 mm x 180 mm (with batt.)
Weight	Approx. 198g (without batt.) Approx. 440g (with batt.)

Specifications may be changed without prior notice.

Authorised Distributor :