




JTMD-201M

DIGITAL THERMOMETER

| | | |
|--|--|--------------------------|
| Measurement: | 32.0°C~43.0°C | 89.6°F~109.4°F |
| Accuracy: | ±0.2°C - 32.0°C~35.4°C | ±0.3°F - 89.6°F~95.7°F |
| | ±0.1°C - 35.5°C~42.0°C | ±0.2°F - 95.9°F~107.6°F |
| | ±0.2°C - 42.1°C~43.0°C | ±0.3°F - 107.8°F~109.4°F |
| | when measured at a standard room temperature of 23°C (73.4°F) in the waterbath | |
| Display unit: | 0.1°C | 0.1°F |
| Storage (transport) condition: -25°C~55°C (-13°F~131°F) /30% - 85% RH | | |
| Operation condition: 10°C~35°C (50°F~95°F) /10% - 95% RH | | |
| Battery: One 1.55V button size battery (LR 41) | | |
| Power consumption: 0.15 MW in measurement mode | | |
| Battery life: More than 200 hours of continuous operation | | |
| Weight: Approx. 20 grams including battery | | |
| Dimensions: 126mm x 22mm x 11.7mm | | |
| Applied part:  type BF | | |
| Degree of protection from liquid and solid particle: IP22 | | |
| The digital thermometer is waterproof | | |



Fast
Accurate



Adult / Child

Authorised Distributor:

