





JTMI-603M NON CONTACT FOREHEAD THERMOMETER



Measuring range:

Human Body: 10°C~50°C (50°F~122°F)

Object : 0°C~100°C (32.0°F~212.0°F)

Measuring accuracy:

Human Body: 36°C~39°C (96.8°F~102.2°F); ±0.2°C
outside of range: ±0.3°C

Object : ≤ 20°C: ± 1°C

> 20°C: ± 5%

Temperature unit: °C and °F selection

Measuring distance: 0.5 ~ 3 cm

Memory: 30 measurements

Operating environment:

16°C~40°C (60.8°F~104°F)

with relative humidity up to 85% (non condensing)

Storage/ Transportation environment:

-20°C~+50°C (-4°F~+122°F)

with relative humidity up to 85% (non condensing)

Display resolution: 0.1°C or °F

Power supply: 2 x 1.5V AAA size alkaline batteries

Weight : approx. 80g (without battery)

Dimensions :

approx. 123.0mm × 50.5mm × 49.0mm (L x W x H)



Adult / Child

Authorised Distributor: